

# Exploratory Data Analysis

Stage 0

# Raw Data

- 25 Hz RT Logs



	HR-PCI	AMI/CGS
Cases	83 (82*)	92 (91**)
Total (min)	251,982	401,780

\*HR-PCI 1991507: empty folder.

\*\*AMI/CGS 2006428: blob not found in csv.


- Case description

- Distribution of RT Log Length



# Variables & Features



- 44 Cols → 18 Variables (Corr)



- 9 Features

Previously	Now
<u>AoP</u> = dP + LVP	<u>Pump Pressure</u> , dP
Motor Current, Motor Speed	Pump Current, Pump Speed, Pump Flow
	Purge Pressure, Purge Flow
	EventSet, EventClr

# Sample RT Logs




# Preprocess

Stage 0

## HR-PCI RT Logs

Patient 1



Patient 2




...

Patient 82

## HR-PCI RT Logs

Patient 1



Patient 2




...


Patient 82

## AMI/CGS RT Logs

Patient 1



Patient 2



...

Patient 91

## Step 1: Concat Consecutive RT Logs

Patient 1



Patient 2



\* Non-consecutive due to missing file  $\Rightarrow$  Treat as 2 logs.

...

Patient 91



## Step 2: Remove Invalid Data

Patient 1



Patient 2




\* Non-consecutive due to removed alarm (event), out of range data ( $[20, 200]$  mmHg)  $\Rightarrow$  Treat as 2 logs.

...

Patient 91

## Step 3: Convert RT Logs to AT Logs

Patient 1




\* Average every 250 data points. Then rolling mean (smoothing) every 5 resulting data points.




## Step 4: Sliding Window Generation of Samples


Patient 1




\* Samples of length 20 minutes.



\* Different forecasting horizons.



\* Train : Validate : Test = 3 : 1 : 1





## Step 5: Trend Classification

\* Different forecasting horizons.






Current MAP (mmHg)	Significant Change
< 60 (critical)	5 mmHg
60 - 80 (managed)	5 mmHg
> 80 (normal)	10 mmHg





# Example Trends

I




D



# Linear Fit Baseline

- Slope-based benchmark analysis to the result previously reported.



	RMSE (mmHg)		RMSE (mmHg)
IDS	9.089	I	12.272
ID	12.614	D	12.728
IS	7.208	S	2.014
DS	7.446		

\* I = increasing, D = decreasing, S = stationary.